## **Quality Mixer**



QA Mixer works against ergonomic injuries in hand wrists and shoulders which happens during manual operation. QA Mixer strips the tube and mix the bag automatically. QA Mixer is adjustable from 3-10 stripping cycles and 3-10 mixing cycles. For example the tube will be stripped 5 times and the bag will be mixed 10 times between each strip. QA Mixer is suitable for all blood bags on the market.

## **Technical data:**

Mains unit inlet 24 V

Mixing cycle 11 cycles per minute

Dimensions 290 (L) x 253 (W) x 150 (H) mm.

Weight 5 kg Max Weight 500 c

Power supply unit inlet 230 VAC + - 10%, outlet 24 V

Symbol CE marked

## **LJUNGBERG & KÖGEL AB**

Box 1032, 251 10 Helsingborg, Sweden Tel: +46 (0) 42-139860 Fax: +46 (0) 42-132181

Fel: +46 (0) 42-139860 Fax: +46 (0) 42-132181 e-mail: info@ljungberg-kogel.com